

# Blood and plasma bag labels



Blood and plasma bags in refrigerators and freezers pose unique labelling challenges, including reliable low-temperature performance and chemical exposure.

More than one million human tissue samples are collected, stored and transported each year – at temperatures as low as -196°C. Our comprehensive low temperature portfolio is designed for bottles, test-tubes, vials and other containers used in hospitals, clinics and research labs. It enables many printing technologies, across applications ranging from blood bags through to cryogenic and laboratory packaging.

Labels comply with FDA and European food regulations, and have been approved by ISEGA testing institute according to DIN ISO 3826 for “Plastic Collapsible Containers for Human Blood and Blood Components”.

## Key features

- ▶ Extreme low cryogenic temperatures: -196°C
- ▶ Steam autoclave sterilisation (up to +121°C)
- ▶ Good chemical resistance
- ▶ UV flexo, thermal transfer, direct thermal and UV inkjet
- ▶ White, opaque and clear versions

## Application areas

- ▶ Blood bag and plasma bag labelling
- ▶ Human and animal health: blood fractioning, stem cells, serum
- ▶ Hospital, clinic and laboratory labelling
- ▶ Deep freeze applications
- ▶ Re-labelling with opaque labels

## Product information

Product Code	Product Description	Service temp. range	Minimum application temp.	Multi-layer construction	Autoclave sterilization	DIN ISO3826 compliance	FDA175.105 compliance	Printability		
								Thermal transfer	Direct thermal	Laser Print

## Primary and secondary application

AD980	TRANSFER PP AL171-BG40WH	-80°C to +140°C	+10°C		✓	✓	✓	✓		
AA134	PE105 MATT WHITE AL171-BG40WH	-80°C to +140°C	+10°C		✓	✓	✓			✓
AA437	MC WETT STRENGTH-AL171-BG40WH	-80°C to +140°C	+10°C		✓	✓	✓	✓		

## Secondary application

AQ665	PB PP NG TOP WHITE 2XS2060NP-BG40BR/BG40WH	-40°C to +70°C	0°C	✓		✓	✓	✓		
AX238	PB PP NG TOP WHITE 2XC2020P-BG40BR/BG40WH	-50°C to +80°C	-20°C	✓		✓	✓	✓		
AX583	THERMAL TOP S8 FSC C2020P-BG40WH	-50°C to +80°C	-20°C			✓	✓		✓	
AO525	TRANSFER PE HD C2020P-BG40WH	-50°C to +80°C	-20°C			✓	✓	✓		
AN038	TRANSFER PEHD S2060NP-BG40BR	-40°C to +70°C	0°C			✓	✓	✓		
AQ596	PP NG TOP WHITE S2060NP-BG45WH	-40°C to +70°C	0°C			✓	✓			

- The product codes highlighted by the green mark are products from the ClearIntent™ portfolio.



Low  
MOQ



Quick  
Quote



Fast  
Delivery

**South Asia Pacific and Sub-Saharan Africa**  
460 Alexandra Road  
PSA Building, #28-02/03  
Singapore 119963  
+65.6349.0333

**North Asia**  
5th Floor, Hongye Park  
1801 Hongmei Road,  
Xuhui District 200233,  
Shanghai, China  
+86 21 33951888

**Europe**  
Willem Einthovenstraat 11  
2342 BH Oegstgeest2300  
AA Leiden  
The Netherlands  
+31 85 000 2000

**Latin America**  
Rodovia Vinhedo-  
Viracopos, KM 77  
CEP 13280-000  
Vinhedo - SP, Brazil  
+55 19 3876-7600

**North America**  
8080 Norton Pkwy  
Mentor, OH 44060  
800.944.8511

DISCLAIMER – All Avery Dennison statements, technical information and recommendations are based on tests believed to be reliable but do not constitute a guarantee or warranty. All Avery Dennison products are sold with the understanding that purchaser has independently determined the suitability of such products for its purposes. All Avery Dennison's products are sold subject to Avery Dennison's general terms and conditions of sale, see <http://terms.averydennison.com>.

© 2018 Avery Dennison Corporation. All rights reserved, Avery Dennison and all other Avery Dennison brands, this publication, its contents and product names and codes are owned by Avery Dennison Corporation. All other brands and product names are trademarks of their respective owners. This publication must not be used, copied or reproduced in whole or in part of purposes other than marketing by Avery Dennison.



Label and  
Packaging Materials

label.sapssa.averydennison.com